

ENDOZIME SPONGE® *Mini*

The Mini Endozime Sponge® is contoured to safely fit around small flexible and rigid scopes and tubular medical instruments with outer diameters of 2 to 6 mm. Providing a quicker and more thorough pre-cleaning.




Cat No 345MINI
Unit of issue - 100

The Mini's unique smaller diameter opening helps protect delicate small scopes from aggressive cleaning where too much pressure may be applied during wiping. Saturated with Ruhof's enzymatic cleaner, Endozime® and designed to start cleaning on contact.

Here are some of the features and benefits of this new sponge:



Easy to use – the Endozime Sponge Mini is specifically contoured to fit 2 mm to 6 mm endoscopes and other tubular instruments.



- Convenient and safe – each Endozime Sponge Mini is packed individually for singular use.
- Uniquely designed to thoroughly clean all surface areas of rigid or flexible endoscopes including flexible cystoscopes, bronchoscopes, pediatric scopes, and ENT scopes etc.
- Hospital-proven enzymatic detergent action – safely and effectively removes organic and proteinaceous materials, significantly reducing micro-organisms.
- Increases effectiveness of manual scope cleaning as speed and simplicity of use allows you to remove gross contaminants while scopes are still wet.
- Does not adversely affect delicate endoscopic parts, such as lenses, lighting fixtures or air/water nozzles.
- 100% biodegradable.
- Eliminates biological odours on contact.