

LIQUIZIME®

MULTI-TIERED ALKALINE ENZYMATIC DETERGENT SPECIFICALLY DESIGNED FOR USE IN WASHER-DISINFECTORS, TUNNEL WASHERS AND CART WASHERS.

Liquizime is a premium detergent that utilizes the power of enzymes in an alkaline solution.

With its multi-tiered formulation of enzymes, Liquizime removes all blood, fat, carbohydrates, starches and protein, while its powerful alkaline solution removes scale and mineral deposits left from hard water. Liquizime's unique sequestering and emulsifying ability assures complete cleanliness and free rinsing of surgical instruments, trays and trolleys.

Description

One-step cleaning and scale prevention.

Liquizime is alkaline (ph 9 to 10) yet does not require neutralizer.

Bio-degradable.

Instantly and constantly active (sequestering).

Free Rinsing, **Liquizime** leaves no residue.

Available in:

Cat No 345LQZ	4 x 4 Litres
Cat No 345LQZ10L	2 x 10 Litres
Cat No 345LQZ20L	1 x 20 Litres